

Designed in Europe Manufactured and assembled in Canada.

The risk of back pain and other musculoskeletal disorders (MSDs) is a very real issue for care givers.

Consider the amount of physical effort that is required and the awkward positions that a care giver may have to get into when transferring a client from a sitting to a standing position.

The use of a **rotunner sit to stand transfer aid** can remove the need for manual handling and lower the risk of injury to the care giver, while providing a safer, more dignified and comfortable transfer solution for the client.



Three Year Warranty

rotunner

THREE YEAR WARRANTY

Manufactured and tested in Canada to the highest standards by **Jaide Products Inc.** the **rotunner Sit to Stand Transfer Aid's** warranty covers defects in materials or workmanship for a period of three years from the date of original purchase.

During the Warranty period **Jaide Products Inc.** will replace or repair, at no cost, any part of the **rotunner** which proves to be defective due to manufacture defects or poor workmanship. If **Jaide Products Inc.** is unable to provide a replacement part or undertake a suitable repair or if such repair is impracticable or cannot be made in a timely fashion, **Jaide Products Inc.** reserves the right to replace the whole unit.

This Warranty does not cover:
Wear and tear, improper usage, misuse or abuse, negligence, accident or alteration.



32 Seapark Dr, St. Catharines,
Ontario, CANADA L2M 6S6



Tel: 289-968-5968

Email: info@jaideproducts.com

rotunner



MADE IN CANADA

roturner - A Dignified Sit to Stand Transfer Solution

Easy • Safe • Comfortable

The **roturner's** sturdy frame instils confidence in the client as they make the transition between a bed, seat, commode or wheelchair. The wide turntable provides a stable and safe area that the client can step onto while the transfer takes place. It rotates between five different positions with a locking device to ensure client safety.


The client's risk and fear of falling is diminished with the curved support frame providing added security and peace of mind. The client is actively involved in the transfer, which will improve their circulation so they can enjoy a sense of independence rather than just being passively lifted from one place to another.

Standard frame roturner

Overall height - 40"/1020mm
Handle height - 30.5"/775mm
Width between handles - 17"/432mm
Overall width (at base) - 26"/660mm
Weight capacity - 551lbs /250kg

Bariatric frame roturner

Overall height - 40"/1020mm
Handle height - 31"/780mm
Width between handles - 20.5"/520mm
Overall width (at base) - 26"/660mm
Weight capacity - 551lbs /250kg



There really is no maintenance required except for cleaning. The **roturner** is powder coated which makes cleaning easy. Simply hand clean with a non- abrasive cloth using mild detergents or disinfectants.